

The average age was 62.3 years in the Supartz®/Supartz FX® group and 74% of the subjects were female. The average age was 61.9 years in the GenVisc® 850 group and 80% of the subjects were female. In both groups, the average weight was 66 kg. There were no statistically significant differences in demographic characteristics.

Primary Effectiveness: In the full analysis set (FAS) population, the VAS pain on movement of the Supartz®/Supartz FX® group at week 6 decreased by 48.0±23.39 mm compared to baseline, and that of the GenVisc® 850 group decreased by 49.2±21.50 mm. The difference between the two groups was not statistically or clinically significant (P>0.05). These analyses are shown in Table 4.

Table 4: VAS pain on movement (mm) and baseline variations (Week 6 – Baseline) (FAS)

	Supartz®/ Supartz FX® group	GenVisc® 850 group	Statistics	P value	Method
N	113	116	0.403	0.688	t test
Mean±SD	−48.0±23.39	−49.2±21.50			
95%CI(Lower-Upper)	−52.33—43.61	−53.12—45.21			
Min-Max	−95.0-1.00	−90.0-17.00			
Median	−50.00	−50.25			

Meta-Analysis of 4 Studies Using Bayesian Modeling

To further support the clinical similarities between GenVisc® 850 and Supartz®/Supartz FX®, a prospective meta-analysis of the pivotal studies for both products was undertaken using Bayesian longitudinal modeling. The studies analyzed include:

- For GenVisc® 850, two saline-controlled studies (AMELIA¹ and Blanco²); and
- For Supartz®/Supartz FX®, two saline-controlled studies conducted in Australia³ and Sweden⁴.

Primary Objectives of the Bayesian Analysis:

- Supartz®/Supartz FX® is superior to PBS. The null hypothesis is that PBS is superior to Supartz®/Supartz FX®. Rejection of the null hypothesis will in effect validate the statistical approach and modeling as it duplicates the results of the approved PMA for Supartz®/Supartz FX®.
- GenVisc® 850 is superior to PBS. The non-inferiority margin for addressing this objective is 4 mm.

Supporting Objective of the Bayesian Analysis:

GenVisc® 850 advantage over PBS is non-inferior to Supartz®/Supartz FX®'s advantage over PBS. The non-inferiority margin for addressing this objective is 4 mm.

Results of Meta-Analysis of 4 Studies Using Bayesian Modeling

Primary Analysis: For the primary analysis, which pools all data from post-baseline visits for all treatments in all studies, the estimated between-study variability (τ) was examined and found to be acceptable for superiority and non-inferiority assessments.

The Gelman-Rubin convergence statistic was very close to 1, thus indicating convergence of the sampler. Overall the model fits the data well.

Primary Analysis

The posterior probability of superiority of GenVisc® 850 vs. PBS is 79% at week 30 (mean 6.88 mm advantage) thus giving confidence that GenVisc® 850 is superior to PBS up to 30 weeks.

Secondary Analysis

For the primary and supporting analyses, differences in mean change from baseline between GenVisc® 850 and PBS were examined. Paucity of data towards the end of the time interval causes an increase in variance and therefore the posterior probability of non-inferiority does not increase. The posterior mean difference between effect GenVisc® 850 and Supartz®/Supartz FX® was always below the non-inferiority margin with a posterior probability of 50%, but the scarcity of data limits the ability to declare non-inferiority of GenVisc® 850 to Supartz®/Supartz FX® for the interval extending to 30 weeks.

Further details of the primary and secondary analysis assessments are provided below in Table 5 and Figure 1.

Table 5. Posterior Probabilities for Main Analysis of GenVisc® 850 Superiority

Objective	π	Probability
GenVisc® Superiority	π ₁	79

The results of the longitudinal analyses are presented in support of the observation that GenVisc® 850 is superior to PBS across time. There is a good linear fit of the data to the model demonstrating increasing mean differences between GenVisc® 850 and PBS through 30 weeks, Figure 1.




Figure 1. Treatment Difference Estimates Across Time Using a Linear Trend Longitudinal Model.

DETAILED DEVICE DESCRIPTION

Each 3mL prefilled syringe of GenVisc® 850 contains:

Sodium Hyaluronate	25.0mg
Sodium Chloride	21.3mg
Disodium Phosphate Dodecahydrate	1.5mg
Sodium Hydroxide	q.s. to adjust pH
Hydrochloric acid	q.s. to adjust pH
Water for Injection	q.s. 2.5mL

HOW SUPPLIED

GenVisc® 850 is supplied as a sterile, non-pyrogenic solution in 3mL pre-filled syringe.

DIRECTIONS FOR USE

GenVisc® 850 is administered by intra-articular injection. A treatment cycle consists of five injections given at weekly intervals. Some patients may experience benefit with three injections given at weekly intervals. Injection of subcutaneous lidocaine or similar local anesthetic may be recommended prior to injection of GenVisc® 850.

Warning: Do not concomitantly use disinfectants containing quaternary ammonium salts for skin preparation because sodium hyaluronate can precipitate in their presence.

Precaution: Do not use GenVisc® 850 if the package is opened or damaged. Store in the original packaging (protected from light) below 86 °F (30 °C). DO NOT FREEZE. Do not use after expiration date indicated on package. The shelf life is 36 months.

Precaution: Strict aseptic administration technique must be followed.

Precaution: Remove joint effusion, if present, before injection GenVisc® 850.

To ensure a tight seal and prevent leakage during administration, firmly HOLD the luer lock when removing the tip cap and attaching the needle. Take care not to rotate the luer lock which can lead to loosening of the hub. Inject GenVisc® 850 into the knee joint through a 21-23 gauge needle.

Inject the full 2.5mL in one knee only. If treatment is bilateral, a separate syringe should be used for each knee.

Precaution: The prefilled syringe is intended for single use. The content of the syringe must be used immediately once the container has been opened. Discard any unused GenVisc® 850.

Please contact your local representative, Customer Care, or scan the QR Code below.



Instructions for Use and Symbol Glossary
eifu.avanos.com

Avanos Medical Sales, LLC,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.
In USA, 1-844-428-2667. www.avanos.com

Meiji Pharma Spain, S.A.
Avda. de Madrid, 94,
28802 Alcalá de Henares,
Madrid, Spain.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. GenVisc® 850 is sold outside the U.S. under the branded name Adant®, a registered trademark of Meiji Seika Pharma Co., Ltd. Supartz and Supartz FX are Registered Trademarks of Seikagaku Corporation. Tylenol is a Registered Trademark of Kenvue Inc. © 2018 AVNS. All rights reserved. 2024-01-10 15-M1-1972-01



Drug therapy

- Pain relievers such as acetaminophen and narcotics
- Drugs that reduce inflammation and pain (signs of inflammation are swelling, pain or redness), such as aspirin and other nonsteroidal anti-inflammatory drugs (NSAIDs) such as ibuprofen, naproxen and Celebrex
- Steroids injected directly into your knee

THINGS YOU SHOULD KNOW ABOUT GENVISC® 850.

- GenVisc® 850 is only for injection into the knee, performed by a doctor or other qualified health care professional.
- If you have pain in both knees, the effects of pain relief on the knee opposite to the one injected have not been evaluated. Your doctor may recommend that both knees be injected.
- After you receive the injection, you may need to avoid strenuous activities such as jogging, tennis, heavy lifting, or standing for a long time for approximately 48 hours.
- If any of the above symptoms or signs appear after you are given GenVisc® 850, or if you have any other problems, you should call your doctor.

WHY IS MY DOCTOR RECOMMENDING I RECEIVE GENVISC® 850?

GenVisc® 850 is indicated for patients with osteoarthritis knee pain who do not obtain adequate relief from simple painkillers, like acetaminophen (Tylenol®), or from exercise and physical therapy. Ask your doctor for additional information and discuss your treatment options.

IS GENVISC® 850 APPROVED FOR USE IN THE U.S.?

GenVisc® 850 was approved for use in the United States in September 2015. GenVisc® 850 has demonstrated a very favorable safety profile, with over 35 million doses administered worldwide since its first approval, and has been studied in over 30 clinical trials for use in relieving osteoarthritis knee pain.

WHAT IS GENVISC® 850 MADE FROM?

Hyaluronic acid is a natural chemical found in almost all species and various parts of your body. It is in high amounts particularly in joint tissues and in the fluid that fills the knee joint space (synovial fluid). GenVisc® 850 is made from an extraction and purification of hyaluronic acid from fermentation of bacteria that make hyaluronic acid identical in chemical composition to human hyaluronic acid.

WILL MY INSURANCE COVER GENVISC® 850?

Most insurance carriers and Medicare cover GenVisc® 850. The process of obtaining reimbursement varies from plan to plan so talk with your doctor's office and your insurance provider before you begin treatment to find out if GenVisc® 850 is covered.

TREATMENT

ARE FIVE INJECTIONS REQUIRED?

Clinical investigations with GenVisc® 850 have demonstrated that five injections one week apart provide optimal pain relief. Some patients have derived benefit following the third injection. The effectiveness of less than three injections has not been evaluated. You should discuss this with your doctor. Completion of the full injection treatment course is recommended to achieve the greatest therapeutic benefit.

HOW DO I KNOW IF GENVISC® 850 IS RIGHT FOR ME?

A doctor is the best person to advise you on any course of treatment. GenVisc® 850 is approved for the treatment of pain due to osteoarthritis of the knee in patients who do not get adequate relief from simple painkillers or from exercise and physical therapy.

AM I TOO OLD FOR GENVISC® 850 INJECTIONS?

There are no specific precautions or contraindications regarding the use of GenVisc® 850 in older patients.

CAN I TAKE OTHER MEDICATIONS WITH GENVISC® 850?

Because GenVisc® 850 is injected directly into the joint rather than administered systemically and there are no known drug interactions, you may be able to take other medications with it. It is not known to interfere with other pain relievers and anti-inflammatory drugs. The safety and effectiveness of the use of GenVisc® 850 at the same time as other intra-articular injections, such as steroids, have not been established. Because it is administered locally it is unlikely that you will observe any benefit in joints not specifically injected. You should discuss all of your current medications and vitamin supplements with your physician on a periodic basis.

CAN I RECEIVE GENVISC® 850 IN BOTH OF MY KNEES?

Yes, if both knees have pain. GenVisc® 850 treatment may be given in both knees simultaneously or separately, according to your physician's recommendations. Because GenVisc® 850 is only injected locally into the knee joint and not administered systemically, it is not likely that injections into one knee will have an effect on the opposite knee if it is not injected. In addition, it is unknown whether a good response in one knee is a good predictor of an equivalent response in the opposite knee should you decide to have this therapy.

SHOULD I MODIFY MY LEVEL OF PHYSICAL ACTIVITY AFTER RECEIVING GENVISC® 850 TREATMENT?

It is recommended that you avoid strenuous activities such as jogging, tennis, and heavy lifting for at least 48 hours after receiving an injection. Your doctor will advise you about the level of activity that is right for you, but in general, patients are able to maintain their normal daily activities after receiving GenVisc® 850 treatment.

CAN GENVISC® 850 BE USED IN JOINTS OTHER THAN THE KNEE?

The U.S. Food and Drug Administration has only approved GenVisc® 850 for use in the treatment of OA of the knee. The FDA has not approved or made an evaluation regarding the safety and effectiveness in other joints.

IS GENVISC® 850 A CURE FOR OSTEOARTHRITIS OF THE KNEE?

There is no cure for osteoarthritis. GenVisc® 850 is a treatment for the symptoms of knee pain of osteoarthritis.

BENEFITS

WHEN CAN I EXPECT TO EXPERIENCE PAIN RELIEF?

Each patient's response to GenVisc® 850 may vary, depending on severity of your OA, degree of pain, and pre-existing medical conditions. In some patients, successful treatment may reduce pain within the first week after treatment begins. However, based upon clinical trials, most patients experienced pain relief after their third injection of GenVisc® 850.

IS GENVISC® 850 TREATMENT EFFECTIVE IN KNEES WITH ADVANCED OSTEOARTHRITIS AND LOSS OF CARTILAGE?

Most of the clinical studies with GenVisc® 850 have selected patients with mild to moderate OA.

WHAT ARE THE BENEFITS OF RECEIVING GENVISC® 850 TREATMENT?

Successful treatment with GenVisc® 850 should reduce pain in an osteoarthritic knee. Because it is a local treatment, GenVisc® 850 should not interfere with any medicine that the patient may take.

HOW LONG CAN I EXPECT THE BENEFIT OF GENVISC® 850 TO LAST?

Each patient reacts differently to GenVisc® 850 treatment. Five injections given at weekly intervals can provide most patients up to 30 weeks of pain relief. The duration of pain relief you experience may vary.

SAFETY

HOW SAFE IS GENVISC® 850?

Extensive safety and toxicity tests were performed on GenVisc® 850. GenVisc® 850 is approved in over 60 countries worldwide. GenVisc® 850 has been in use since 1995 in Japan, 1997 in Europe, and 2015 in the U.S., with over 4 million injections given worldwide. In clinical trials, transient local pain or swelling occurred with some patients with injections of GenVisc® 850. Since introduction into the global market, GenVisc® 850 has not been withdrawn in any market because of safety concerns.

For more Information, please ask your doctor or refer to the Full Prescribing Information for GenVisc® 850. Please contact your local representative, Customer Care, or scan the QR Code below.



Instructions for Use and Symbol Glossary
eifu.avanos.com

Avanos Medical Sales, LLC,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.
In USA, 1-844-428-2667. www.avanos.com

Meiji Pharma Spain, S.A.
Avda. de Madrid, 94,
28802 Alcalá de Henares,
Madrid, Spain.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. GenVisc® 850 is sold outside the U.S. under the branded name Adant®, a registered trademark of Meiji Seika Pharma Co., Ltd. Supartz and Supartz FX are Registered Trademarks of Seikagaku Corporation. Tylenol is a Registered Trademark of Kenvue Inc. © 2018 AVNS. All rights reserved. 2024-01-10 15-M1-1972-01

